

Dux Dental Material Safety Data Sheet

Van R, Cadco, Clive Craig
600 E. Hueneme Road
Oxnard, CA, USA

Phone: 805.488.1122
Fax: 805.488.2266
www.duxdental.com

Product Gingibraid, Gingiknit, Gingiyarn, Unibraid, & Uniknit Alum

Other Name Aluminum Potassium Sulfate Gingival Retraction Materials

REF 13076, 13077, 13078, 13210, 13211, 13212, 13213, 13360, 13361, 13362, 13363, 13365, 13366, 13367, 13368, 13508, 13509, 13510, 13511, 13512, 13513, 13608, 13609, 13610, 13611, 13612, 13613

MSDS Number MS113

Hazardous Ingredients

Principal Hazardous Components	pel / tlv	CAS Number	Percentage
None	NA	NA	NA

Physical and Chemical Characteristics

Appearance	Small cords of varying diameters and constructions impregnated w/ Alum.
Boiling Point	NA
Melting Point	NA
Vapor Pressure	NA
Evaporation Rate	NA
Odor	None
Vapor Density	NA
Specific Gravity	NA
Water Solubility	Yes (Alum)

Fire & Explosion Data

Explosion	Not considered to be an explosion hazard.
Fire	Not considered to be a fire hazard.
Fire Ext. Type	A
Flashpoint	NA
Special Fire Procedures	None
Unusual Hazards	None Known

Reactivity Data

Stability	Stable
Incompatibilities	Sl. corrosive to metals.
Hzd Decomp Products	Oxides of Sulfur
Hazardous Polymerization	Will not occur.

Dux Dental Material Safety Data Sheet

Van R, Cadco, Clive Craig
600 E. Hueneme Road
Oxnard, CA, USA

Phone: 805.488.1122
Fax: 805.488.2266
www.duxdental.com

Product Gingibraid, Gingiknit, Gingiyarn, Unibraid, & Uniknit Alum

Other Name Aluminum Potassium Sulfate Gingival Retraction Materials

REF 13076, 13077, 13078, 13210, 13211, 13212, 13213, 13360, 13361, 13362, 13363, 13365, 13366, 13367, 13368, 13508, 13509, 13510, 13511, 13512, 13513, 13608, 13609, 13610, 13611, 13612, 13613

MSDS Number MS113

Health Hazard Information

Ingestion NA

Inhalation Dust may cause irritation.

Skin/Eye Contact Dust may cause irritation. Flush area w/ water if exposed.

Toxicity Data May cause reversible or irreversible changes to exposed tissue, not permanent injury or death; can cause considerable discomfort. Use in accordance with manufacturers instructions and sound professional

Acute Hazards None

Chronic Hazards None

Carcinogenicity NA, no product ingredient listed by NTP, IARC, or OSHA.

Precautions for Safe Handling and Use

Cleanup NA

Disposal Dispose of in accordance with local, state, and federal regulations.

Precautions As part of good industrial, personal hygiene, and safety practice avoid all unnecessary exposure to the substance and ensure prompt removal from skin, eyes, and clothing.

Occupational Control Measures

Eye Protection NA

Gloves NA

Dust/Particulate Mask NA

Other NA

Additional Information

DOT Hazard Class Not Regulated

DOT Small Qty Exception NA

DOT Other Exception NA

DOT Exception Permitted NA

Allergenic Ingredients This product does not contain Latex, Peanut Products, Gluten, or Red Dye #3.

The information contained herein is provided in good faith and is believed to be correct as of the date hereof. However, the company makes no representation as to the comprehensiveness or accuracy of the information. It is expected that individuals receiving the information will exercise their independent judgement in determining its appropriateness for a particular purpose. Accordingly, we will not be responsible for damages of any kind resulting from the use of or reliance upon such information. No representations or warranties, either express or implied, of merchantability, fitness for a particular purpose or of any other nature are made hereunder with respect to the information set forth herein or to the product to which the information refers.